



Manufactured in U.S.A. by:

# SCI-PHARM



Sci-Pharm Catalog No.  
70-125 B

## GEL ETCH™

### TECHNICAL BULLETIN

#### Instructions


# GEL ETCH™

## ECONOMY-SIZE DISPENSING SYRINGE

Cat. No. 70-125 B

### KIT CONTAINS:

- 60g Etching Agent
- 10 Applicator Syringes
- 60 Disposable Dispensing Tips



#### CAUTION

Contains Phosphoric Acid

### OUTSTANDING FEATURES OF THE MATERIAL

Scientific Pharmaceuticals' *Gel Etch*™ is a sophisticated dental device for the conditioning of tooth surfaces prior to restorative and preventive procedures, including cementing, bonding, sealing, and filling. The use of special additives assures good wetting of tooth enamel and dentin and a stable consistency that resists drying. The FD&C colorant provides good contrast with mucosa and allows for better control of etched areas, while allowing for easy rinsing of the material from teeth and soft tissues.

Packaged in our economy-sized dispensing syringe, it provides a convenient and economical solution to providing an advanced etchant formulation for busy group practices, multichair operations, health institutions, charitable and special programs, dental schools, etc. Small applicator syringes are easy to fill from the large dispensing syringe and, in combination with disposable tips, provide a safe, precise, and extremely convenient delivery system.

#### Instructions for Filling Applicator Syringes from the 60g Dispensing Syringe

Remove caps from both syringes; attach the small syringe to the adapter of the large syringe with the piston of the small syringe pushed all the way in; ensure that the small syringe is well-locked onto the large syringe by twisting the small syringe to the right; slowly exert pressure on the piston of the large syringe; material should flow smoothly from the large syringe into the small syringe - if resistance to flow is encountered, pull the piston of the small syringe simultaneously; remove small syringe; clean any excess material with paper towels. If contact with hands occurs, wash with water immediately; replace caps on both syringes.

#### CAUTION

Federal law restricts this device to sale  
by or on the order of a dentist.

